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Multimat[®]2 Touch

Multimat[®]2 Touch+Press

Instructions for use
Instructions de service
Istruzioni per l'uso
Manual de instrucciones

Multimat®2 Touch

Multimat®2 Touch+Press

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Stand der Information 27.07.04

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Software V1.01 and higher

Information current as of July 27, 2004

All parts of this operating manual are copyrighted.

Copying or circulating this text in whole or in part is only legal with written permission of the DeguDent GmbH.

Violation obliges to compensation. Multimat® is a registered trade mark of the DeguDent GmbH.

Subject to technical and design changes.

© **DeguDent GmbH, Hanau-Wolfgang**

Introduction

Dear customer,

Thank you for purchasing the Dentsply **Multimat®2 Touch**.

This furnace represents a further development of the **Multimat®2 Touch** series. Due to its technical capabilities, it is excellently qualified to match future developments in dental ceramics. By choosing from various firing modes, this furnace allows you to process materials of different manufacturers largely in familiar way. The simple, menu-driven operation of the **Multimat®2 Touch** keeps training time at a minimum.

The **Multimat®2 Touch** is equipped with 60 fixed and special programs as well as an additional 300 freely programmable programs. The touch-screen allows for direct and quick selection of all functions. All firing data are shown on a graphic color display. A 32-bit microcontroller together with a real-time operating system controls the measuring and regulating processes. The Touch+Press version allows pressable dental ceramics to be processed in addition to the usual firing options. The instrument complies with all valid EU directives and VDE/UL safety regulations.

Please read the operating manual prior to start-up! It contains important information on safety, operation and maintenance, to protect you and prevent damage to the instrument.

Please keep this operating manual and include it when passing this instrument on to another owner.

Multimat®2 Touch

Contents

INTRODUCTION	38
WARNING NOTES	42
 1 Parts list	 43
1.1 PARTS LIST	43
1.2 TECHNICAL DATA	44
1.3 ENVIRONMENT CONDITIONS	44
 2 Safe use	 45
2.1 USE WITHIN SPECIFICATIONS	45
2.2 HAZARDS AND SAFETY NOTES	45
 3 Setup and first operation	 45
3.1 UNPACKING	45
3.2 CHECK ACCESSORIES	45
3.3 SETUP	46
3.4 INITIAL OPERATION	46
3.4.1 Start screen display	46
3.4.2 Language selection	46
3.4.3 Introduction	46
3.4.4 Installation and safety notes	46
3.4.5 Functional test	46
3.4.6 Basic settings	47
3.4.6.1 Temperature mode	47
3.4.6.2 Night/Standby temperature	47
3.4.6.3 Vacuum mode	47
3.4.6.4 System time	47
3.4.6.5 Date	47
3.4.6.6 Acoustic signal	47
3.4.6.7 Idle temperature	47
3.4.6.8 Upper temperature limit	47
3.4.6.9 Data output	47
3.4.6.10 Screen brightness	48
3.4.7 Ceramic type	48
3.4.8 Program: special programs	48
3.5 TEST PROGRAMS	48
3.5.1 Test program Dentsply / Ceramco / External	48
3.5.2 Start test program	48

4	Practical application: introduction	49
4.1	FIRING PARAMETER LIMIT VALUES	49
4.2	SCREEN	50
4.3	SCREEN FUNCTIONS	50
4.4	SENSOR KEYPAD FUNCTIONS	50
4.5	MAIN MENU	51
4.5.1	Configuration	51
4.5.2	Furnace parameters	52
4.5.2.1	Reset calibration offset	52
4.5.2.2	Vacuum test	52
4.5.2.3	Reset language	52
4.5.2.4	Heating % (power factor)	52
4.5.2.5	Reset hours of operation (muffle)	52
4.5.2.6	Statistical data	52
4.5.3	Smart <i>Media</i> Card	52
4.5.4	Delete programs	52
4.5.4.1	Single program	52
4.5.4.2	All individual programs	53
4.5.4.3	Delete fixed and special programs	53
5	Create individual programs	53
5.1	FIRING METHOD	53
5.2	PROGRAM TYPE	53
5.3	PROGRAM LIST	53
5.4	TEXT INPUT	53
5.5	FIRING PARAMETERS	53
5.5.1	Set preheating temperature	54
5.5.2	Set predrying time	54
5.5.3	Set drying time	54
5.5.4	Set preheating time	54
5.5.5	Set vacuum level	54
5.5.6	Set heating rate	54
5.5.7	Set firing temperature	54
5.5.8	Set vacuum	54
5.5.9	Set firing time	54
5.5.10	Set tempering temperature	54
5.5.11	Set tempering time	54
5.5.12	Set cooling stage	54
5.6	SAVE PROGRAM	54
5.7	START PROGRAM	55
5.8	STANDARD FUNCTIONS	55
5.8.1	Opening a program	55
5.8.2	Saving a program	55
5.8.3	Changing a program	55
5.8.4	Changing a program during the firing sequence	55
5.8.5	Copying a program	55
5.8.6	Deleting programs	56
5.9	PRESSING	56

6	Special functions	56
6.1	NIGHT/STANDBY	56
6.2	QUICK COOLING	56
6.3	VIEW FUNCTION	57
6.4	POWER FAILURE BACKUP	57
6.5	TESTING THE FIRING CHAMBER TEMPERATURE	57
6.6	NIGHT MODE	58
6.7	SELECT ACTIVATION TIME	58
6.8	VACUUM TEST	58
6.9	PRINT/PC (PC PROGRAM OPTIONAL)	58
6.10	ACOUSTIC SIGNALS	58
6.11	SOFTWARE UPDATE	58
6.12	PROGRAM SORTING BY NAME OR PROGRAM NUMBER	59
6.13	CHANGE STANDARD DRYING POSITION AND STEPS	59
6.14	CHANGE TEMPERING POSITIONS	59
6.15	DISPLAY DATE OR TIME	59
6.16	SMART MEDIA CARD	60
6.16.1	Program transfer furnace → Smart Media Card	61
6.16.2	Program transfer Smart Media Card → furnace	61
7	Service and maintenance	61
7.1	LIFT	61
7.2	FIRING SOCKET/PRESSING SOCKET AND FIRING SOCKET CARRIER	61
7.3	VACUUM PUMP	61
7.4	CHANGING THE FIRING MUFFLE (TRAINED ELECTRICIANS OR SERVICE PERSONNEL ONLY!)	62
7.5	REPLACING THE CONTROL UNIT (TRAINED ELECTRICIANS OR SERVICE PERSONNEL ONLY!)	62
7.6	INSPECTIONS	62
7.7	CLEANING NOTES	62
8	Faults and remedies	63
9	Preset programs	65
9.1	FIRING TABLE FOR CERAMCO 3 METAL CERAMICS (DENTSPLY MODE)	65
9.2	FIRING TABLE FOR FINESSE METAL CERAMIC (DENTSPLY DETREY MODE)	66
9.3	9.3 FIRING TABLE FOR FINESSE ALL CERAMIC (DENTSPLY MODE)	67
9.4	FIRING TABLE FOR FINESSE, FAC, CERAMCO II (DENTSPLY CERAMCO MODE)	67
10	Warranty	68
11	Declaration of Conformity	69
12	Dentsply Subsidiaries	172

Figures

1	Parts list	43
2	Muffel	43
3	Keypad	43
4	Language selection menu	46
5	Introduction screen	46
6	Installation menu	46
7	Functional test menu	46
8	Basic settings menu	47
9	Ceramics types menu	48
10	Special programs menu	48
11	Test program	48
12	Keypad	50
13	Main menu	51
14	Configuration menu	51
15	Furnace parameters menu	52
16	Program deletion menu	52
17	Firing modes menu	53
18	Program types menu	53
19	Metal ceramics menu	53
20	Entering characters	53
21	Form for entering firing data	53
22	Pressing program menu	56
23	+3.3 V Smart <i>Media</i> Card with 128 MB of memory	60
24	+3.3 V Smart <i>Media</i> Card with write protection in place	60
25	Inserting the Smart <i>Media</i> Card	60
26	Smart <i>Media</i> Card eject button	60
27	Smart <i>Media</i> Card menu	60
28	Smart Card, front and rear view	61

Warning notes

Hazard warning:

This symbol indicates particular hazards, where non-compliance with the instructions can cause injury and accidents.

Hot surface warning:

This symbol indicates a hot surface on the instrument.


***Texts in italics only apply to the furnace version
Multimat2 Touch+Press.***

1 Parts List


1.1 Parts list


Figure 1: Parts list

- 01. Compression hood
- 02. Cover
- 03. Cooling jacket
- 04. Support
- 05. Firing platform
- 06. Firing platform holder
- 07. Placement area
- 08. Control housing
- 09. Smart *Media* Card slot
- 10. USB slave interface (PC interface)
- 11. RS 232 interface (printer/PC interface)
- 12. RJ 45 interface (PC network)/RJ12 (for ODS)
- 13. Touch screen
- 14. Control
- 15. Fuses
- 16. Main power socket
- 17. Vacuum connection
- 18. Vacuum pump socket
- 19. Compressed air connection
- 20. Filter control (no picture)
- 21. Muffle
- 22. Compression socket (no picture)

View of right-hand side

View of left-hand side


Figure 3: Keypad


Figure 2: Muffle

1.2 Technical data

Note:

The instruments Multimat2 Touch und Multimat2 Touch+Press are delivered in one of for 2 different voltage systems, including a power lead (as required for the country in which it is operated):

- 100 – 125 Volt ~ / 50/60 Hz
- 230 – 240 Volt ~ / 50/60 Hz

Check the voltage rating of your instrument on the type plate, which is located on the rear of the instrument. Make sure that the voltage rating and the supplied power cord are appropriate for your local power network.

	Multimat2 Touch	Multimat2 Touch+Press
Power consumption	1,350 W (without vacuum pump)	1,350 W (without vacuum pump)
Pump connection	As mains voltage, max. 2.5 A	As mains voltage, max. 2.5 A
Overvoltage category	II	II
Protection class	I	I
Instrument fuses	2x 16 A medium slow, 6.3 mm x 32 mm, 250 V	2x 16 A medium slow, 6.3 mm x 32 mm, 250 V
Degree of contamination	2	2
Weight	22 kg	25 kg
Dimensions H/W/D	441/320/425 mm closed	593/320/425 mm closed
Dimensions H/W/D	585/320/425 mm open	748/320/425 mm open
Pressing pressure	–	2.7 bar

1.3 Environment conditions

Environment temperature: 2 °C to 40 °C
Relative humidity: 80 % at 31 °C
Altitude: max. 3500 m above sea level

2 Safe use

2.1 Use within specifications

Caution!

The **Multimat2 Touch** was exclusively designed for firing and/or pressing of dental ceramics and is for these uses exclusively. The user is liable for all damages resulting from non-specified use. At temperatures above 1000 °C the working life of the heating muffle and the compression cylinder and the compression valve is reduced. In this case, the aforementioned parts are excluded from warranty. Such use within specification furthermore includes the compliance with all notes in this operating manual as well as all notes in the separate operating manuals of the vacuum pump, the **Online Diagnosis System** and the **Touch+Save PC** documentation.

Repairs and maintenance may only be carried out by Dentsply Service or authorized Dentsply service points.

Operate the touch screen only within the designated sensitive surface. Never use hard or sharp objects to operate the touch screen.

Use the compression platform only for pressing.

2.2 Hazards and safety notes

Caution!

In order to ensure safe operation, please observe the following notes:

- Do not set up furnace and vacuum pump next to heat sources.
- The power consumption of the vacuum pump should not exceed 2.5 A. The powers of muffle and pump add up and can cause overload to the instrument fuses.
- The distance to the next wall should be at least 25 to 30 cm.
- The surface should be non-flammable, and there should be no flammable objects nearby.
- Place the vacuum pump in a well ventilated place. Vacuum pumps lubricated with oil should always be set up in a lower place than the furnace. The back-pressure valve in the tube, however, should be placed higher than the pump.
- Do not touch parts that become hot, particularly the cover.
- Prior to powering up, ensure that the mains power does not exceed the voltage specified on the instrument.
- For 220/230V, use the grounded power lead H05VV-F 3G1,0 with grounded power socket and Schuko-type (grounded) plug St.30 D.
- For voltages 100 – 127 V, use the power lead SJT 3x18AWG 105 °C with UL plug NEMA 5-15P and grounded socket.
- *Furnaces with pressing function must always have connection to compressed air, even during normal vacuum firing. Always maintain the air pressure of 2.7 bar shown on the warning label next to the compressed air inlet.*

- On first-time operation, after any extended rest period as well as at high humidity or low temperatures, it may no longer be possible to generate a sufficient vacuum. In these cases, start program 376.
- If the furnace is under vacuum in “Off” state for an extended period of time, the O-ring of the lift plate may adhere slightly.
- When heating the firing muffle, there may be vibration noises from the heater winding.
- In case of undervoltage, the temperature increase may be delayed.
- Use only original service parts.

Warning!

- The instrument may only be connected to a separate Schuko-type socket which is secured with a slow 16 A fuse. For increased safety, we recommend connecting the instrument to a circuit with separate protective circuit breaker with a release current of 30 mA.
- Always disconnect the instrument from main power when carrying out repairs and maintenance work.
- Never connect the instrument to extension cables or distribution boxes, as there is fire hazard with these devices.
- **Repairs on opened unit under voltage may only be carried out by trained electricians.**
- At least once per year, a specialist must test the protective ground wire.
- In the case of faults or damages that prevent safe operation, the instrument must be secured against unintentional use.
- Do not change the position of the belt tension screw located on the rear column plate.

3 Setup and first operation

3.1 Unpacking

- First, please check the “Shockwatch” label on the carton. If it is colored red, the impact energy during transport was higher than allowed and your unit could be damaged. Ask the transport agent to confirm the triggering of the “Shockwatch” label in writing.
- Open the packaging and unpack the **Multimat2 Touch** carefully.
- Prior to setting up, check the instrument for visible exterior damage. Do not set up or operate a damaged instrument. Please contact us, if the furnace is damaged.

3.2 Check accessories

- Check the delivery on completeness.
Each **Multimat2 Touch** is supplied with:
 - main power lead
 - 1 honeycomb firing goods carrier
 - 1 firing socket
 - 1 pair of tweezers
 - 1 operating manual

Each **Multimat2 Touch+Press** is additionally supplied with:

- 1 compression socket
- 1 filter control with manometer and compressed air tube

For pressing FAC pressable ceramics, you will furthermore need 3.4.2 Language selection aluminum oxide pistons, embedding paper and muffle set.

These parts are included in the FAC complete set. Should you require individual parts, you can order these with the numbers

- D430112 FAC large muffle set
- FAC small muffle set (not available in Germany)
- D430114 FAC aluminum oxide piston
- D430115 FAC embedding paper

- Please retain the transport package.
The supplied package protects the instrument from transport damage. If possible, keep the package material, in case you need to transport the instrument later.

3.3 Setup

- Place the furnace on a suitable surface and allow for sufficient distance to walls (min. 25 cm).
- Place the firing socket on the firing socket holder.
- Connect the power plug of the Dentsply vacuum pump with the vacuum pump socket on the furnace and insert the vacuum tube in the tube connector on the furnace. The arrow on the filter must point into the direction of the vacuum pump.
- If you have a **Multimat2 Touch+Press**, insert the compressed air tube of the filter pressure control into the Schott plug connection of the furnace.
- Connect the filter pressure control with the compressed air feed and set the operating pressure of the pressure furnace to 2.7 bar (the operating pressure of 2.7 bar is set by the manufacturer).

3.4 Initial operation

- Before establishing the power connection, check whether the mains voltage of your power supply corresponds to the mains voltage specified on the type label of your furnace.
- Connect the power lead with the furnace and with a mains socket. The green power LED above the numeric sensor pad illuminates.

3.4.1 Start screen display

The start screen display showing the Dentsply logo, the version number of the software as well as the serial number of the control unit will appear on the screen for 6 seconds.

The firing chamber with a possibly vacuum will be filled with air. The following screen will appear:

3.4.2 Language selection

Figure 4: Language selection menu

Select your language by pressing the appropriate round, yellow sensor button. The sensor button changes color upon pressing. There is no time limit for language selection. The language may also be changed later.

The language selection will be displayed only after resetting the language in the main menu (see section 4.5.2.3) and a complete disconnection from main power.

3.4.3 Introduction

Figure 5: Introduction screen

Press the right arrow to continue to the next screen.

3.4.4 Installation and safety notes

Figure 6: Installation menu

The most important installation notes as well as the note to follow the installation and safety notes in the operating manual are displayed. Confirm with the OK button that you have complied with the installation notes and press the right arrow.

3.4.5 Functional test

Figure 7: Functional test menu

The functional test starts automatically and is documented in the display. The test result is shown in an info field. If there are no errors and all test sections are ticked off, the next screen "Basic settings" is displayed automatically after completion of the test.

If there are errors during the test, the test section is marked with a red X and an appropriate note in the info field is displayed. In this case, confirm the test result with the sensor button “OK”.

The functional test is carried out after any disconnection from main power.

Note: When pressing the left arrow, a safety interval between the screen “Basic settings” and “Functional test” is set to prevent accidental starting of the functional test. You may abort the functional test by immediately pressing the sensor button “OK”.

3.4.6 Basic settings


Figure 8: Basic settings menu

Note:

Several settings of the **Multimat2 Touch** are preset on delivery. You may accept these settings by pressing the sensor button “⇒” or change as described hereafter.

3.4.6.1 Temperature mode

Pressing the yellow sensor button changes the temperature mode between °C and °F.

3.4.6.2 Night/Standby temperature

This function is used to set the temperature of the furnace when switched off but not disconnected from main power. The night/standby temperature prevents moisture to collect in the furnace chamber. This temperature is used in the night mode (see section 4.4) and in standby mode. When switching the furnace off with the “On/Off” switch, the standby function is activated automatically.

The yellow sensor button can be used to set the desired temperature.

1. Press the sensor button (button changes color).
2. Enter the new temperature with the numeric sensor keypad.
3. Temperature >100 °C = night/standby ready
4. Temperature <100 °C = night/standby not ready
5. Press sensor button again (or select next sensor field), and the new temperature value is saved (sensor field changes to yellow again).
6. Wrong entries can be deleted with the sensor button “C”.

3.4.6.3 Vacuum mode

Pressing the yellow sensor button changes the vacuum mode between hPa, inHg and Hg”.

In the modes inHg and Hg”, the pressure is displayed with one decimal.

Note:

hPa → 1013 hPa = standard pressure at sea level
 inHg → 29.9 inHg = standard pressure at sea level
 Hg” → 0 Hg” = standard pressure at sea level (US standard)

3.4.6.4 System time

1. Press the sensor button next to “Time” (button changes color).
2. Enter the 4-digit number for hours/minutes with the numeric sensor keypad.

3.4.6.5 Date

1. Press the sensor button next to “Date”.
2. Enter the 8-digit number for day/month/year with the numeric sensor keypad.

3.4.6.6 Acoustic signal

Press the sensor keypad to switch the acoustic signal on or off.

3.4.6.7 Idle temperature

This function is used to save energy and reduce the surface temperature of the furnace.

Note: The idle temperature is factory-set to 400 °C. It cannot be set to temperatures > 600 °C. The idle temperature must be set to 25 °C lower than the preheating temperature.

1. Press the sensor button next to “Idle temperature”.
2. Enter the new temperature with the numeric keys of the right control pad.

3.4.6.8 Upper temperature limit

This function is used to limit the maximum permissible temperature of the furnace.

Note: It is recommended to limit the upper temperature value, if a heavy stress to the muffle can and shall be prevented.

1. Press the sensor button next to “Upper temperature limit”.
2. Enter the new temperature with the numeric keys of the right control pad.
3. Press the sensor button “Upper limit value” again (or select next sensor field), and the new temperature value is saved (sensor field changes to yellow again).

3.4.6.9 Data output

By pressing the sensor button, the data output changes between “printer”, “PC”, and no output.

If the symbols “printer” or “PC” are activated, all target and actual values are automatically transferred to printer/PC at the end of each firing cycle, as long as a printer or PC is connected.

The theoretical firing data can be displayed at any time by activating the PR/PC sensor button at the bottom of the screen (firing data screen).
If no printer or PC symbol is displayed, the data output is disabled.

3.4.6.10 Screen brightness

The screen brightness can be set with the sensor buttons “+” for increasing the brightness and “-” for decreasing the brightness.

After completion of the basic settings, press the sensor button “⇒”. The furnace changes to the following screen:

3.4.7 Ceramic type


Figure 9: Ceramics type menu

Select the ceramic type from the displayed possibilities by pressing the appropriate sensor button, e.g. special programs.

After selecting the ceramic types “Sinter ceramic”, “Pressable ceramic” and “Special programs”, the furnace changes directly to the appropriate program list.

After selecting the ceramic types “Sinter ceramic”, “Pressable ceramic” and “Special programs”, the furnace changes directly to the appropriate program list.

3.4.8 Program: special programs

Figure 10: Special programs menu

Select the desired program from the program list by pressing the appropriate sensor button, e.g. the test program DeTrey.

The test programs differ only in that they have different entering options for the firing parameters. If you already own a Dentsply furnace, you will quickly recognize the operation of the DeTrey interface. If you already have used Ceramco furnaces, you will probably like the Ceramco test program better. If you work with foreign ceramics, the test program External may be better for you. Give it a try.

3.5 Test programs

3.5.1 Test program Dentsply / Ceramco / External

These test programs will provide you with a first impression of the concrete program functions of the **Multimat2 Touch**. Please see section 6 for details.

3.5.2 Start test program

- Open the firing chamber by pressing the sensor button “↑” on the right sensor pad.
- Select test program.
- Press the sensor button “**Start/Stop**”. The screen program procedure is displayed.


Figure 11: Test program

After starting, the **Multimat2 Touch** automatically changes to the screen program procedure. The status bar above the temperature/time display shows the current program section and the remaining time.

The program starts by heating from the preset idle temperature to the preheating temperature. During this, the firing chamber is open. After reaching the preheating temperature, the program procedure starts with the first program section.

Drying

During this phase, the firing chamber gradually moves downward and the program curve builds up on the screen in dependence of time. Below the program section, the appropriate time is displayed, counting down to zero.

The vertical lift of the firing chamber normally is 150 mm and during compression 158 mm. At 0 mm, the firing chamber is closed. At 150 (158) mm, the firing chamber is in its upper end position. The standard position for drying is at 80 mm. During drying, the firing chamber moves in 9 uniform steps from the upper end position to the standard position for drying. The standard position “Drying 80 mm” and the 9 steps are factory-preset values. The end position and the steps can be changed within these values (see section 6.13 Standard position and Change steps). We recommend keeping these settings, if there are no urgent reasons to do otherwise.

Preheating

The firing chamber moves from the last drying position to the preheating position and the preheating time is displayed counting down to zero.

Time-controlled vacuum (Dentsply firing method)

After completion of the preheating time, the firing chamber closes; the vacuum chamber activates and evacuates the firing chamber until the preset vacuum is reached.

Temperature-controlled vacuum

(Dentsply Ceramco firing method)

Vacuum activation and deactivation is temperature-controlled.

Heating rate

Once the preset vacuum is reached, the temperature increases with the preset heating rate to the firing temperature. Below the program section, the heating time is displayed, counting down to zero.

Vacuum time

After reaching the firing temperature, the firing time under vacuum starts. Below the program section, the vacuum time is displayed, counting down to zero. At the end of the vacuum time, the firing chamber is ventilated.

Firing time

The vacuum time is followed by the firing time at normal pressure, i.e. without vacuum. Below the program section, the firing time is displayed, counting down to zero.

Note: When heating the firing chamber, there may be vibration noises from the heater winding.

Tempering (not programmed in the test program)

The tempering position of 50 mm is a factory-preset value. For other metal ceramics, ask the ceramics manufacturer for the values for tempering position and tempering temperature (see also section 6.14).

Cooling stage (not programmed in the test program)

(See section 5.5.12)

4 Practical application: introduction

4.1 Firing parameter limit values

At the end of the firing time, the firing chamber moves to the upper end position. The end of firing is indicated by a triple acoustic signal and the screen display changes back to the program start interface.


Parameter	Lower limit	Upper limit
Night/standby temperature	101 °C (214 °F)	300 °C (572 °F)
Idle temperature	30 °C (86 °F)	600 °C (1112 °F)
Theoretical value of firing temperature	30 °C (86 °F)	1200 °C (2192 °F)
Actual value of firing temperature	30 °C (86 °F)	1250 °C (2282 °F)
Times	00:00 min or hrs	99:59 min or 17:59 hrs
Controlled heating time	0.1 °C/min (0.18 °F)	120 °C/min (248 °F)
Cooling stages	0	3
Vacuum level	1 hPa (0,1 inHg)	1013 hPa (29,9 inHg)
View	00:01 min	03:00 min
Vacuum ON *	30 °C (86 °F)	1200 °C (2192 °F)
Vacuum OFF *	30 °C (86 °F)	1200 °C (2192 °F)
Firing chamber position		150 mm (<i>compression at 158 mm</i>)
Steps	1	9

* Temperature-controlled vacuum is only selectable in the Ceramco firing method.

Note: Values outside of these limits can neither be set nor started. The entered value is set back automatically to the previous value.

4.2 Screen

Graphic color screen for displaying the firing parameters and the info text. Basic screen layout (English representing all languages):

- | | |
|-----------------------|-------------------------------------|
| 1. Help button | 7. Date/time |
| 2. Program name | 8. Total firing time – theoretical |
| 3. Actual temperature | 9. Firing temperature – theoretical |
| 4. Actual vacuum | 10. Data input |
| 5. Program number | 11. Information field |
| 6. Lift position | 12. Soft-touch sensor buttons |

4.3 Screen functions

1. Help (?)

The help function includes text to provide the user with information in the case of errors. The text included is limited to the most important information.

2. Program name

Screen display for the program name. The program name is taken from the text input field.

3. Temperature – actual

Displays the current temperature in the firing chamber. The temperature unit can be selected in the basic settings section.

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Displays the current underpressure in the vacuum system, from the environment pressure to the preset vacuum level. The vacuum unit can be selected in the basic settings section.

5. Program number

Shows the program number

6. Lift positions

This value indicates the distance of the firing chamber to the firing socket holder.

- If the firing chamber is closed, the value is 0 mm.
- If the firing chamber is opened for **firing**, the value is 150 mm.
- If the firing chamber is opened for **compression**, the value is 158 mm.

7. Date/time

Displays either date or time.

Please select the display mode in the sub-menu “Configuration”.

For documentation printing purposes it is recommended to display the date.

8. Total firing time – theoretical

Displays the estimated addition of all firing times as theoretical value.

9. Firing temperature – theoretical

Displays the current theoretical value for the temperature. The temperature unit can be selected in the basic settings section.

10. Data input

Enter your parameters into these fields. Press the yellow round sensor buttons to start and end the data input.

11. Information field

This screen area displays information during furnace operation.

12. Soft-touch sensor buttons

Press buttons for activating context functions.

4.4 Sensor keypad functions


Figure 12: Keypad

Lift up Power LED Recall keypad

Power LED Lights up, when the **Multimat2 Touch** is connected to main power.

Numeric keypad

For entering numerical values.

Sensor button Recall “R” (start program)

Press the button, enter the program number and press the button again to display the desired program.

Sensor button Lift up (“↑”)

Press the sensor button “↑” to move the firing chamber upwards or stop a downwards movement.

Sensor button Lift down (“↓”)

Press the sensor button “↓” to move the firing chamber downwards or stop an upwards movement.

Sensor button Menu (“M”)

Press this sensor button to display the main menu.

Sensor button Save (“S”)

Is used to save programs (for details see section 5.6).

Sensor button “On/Off”

Press this sensor button to switch the furnace control on or off. The firing chamber needs to be closed to switch the furnace off.

Sensor button “Night Mode”

This function activates the automatic power down of the furnace after firing. The function can be added to the running program at any time by pressing the button with the bed symbol. The bed symbol for Night Mode is displayed in the firing data view and in the progress view.

After the firing is completed, the unit switches off (screen off, muffle off) and the firing chamber remains opened until the Night/Standby temperature is reached. After reaching the Night/Standby temperature, the firing chamber closes. Pressing the On/Off button switches the furnace on again.

Sensor button Quick cooling “☒”

Press the sensor button “☒” to manually activate the quick-cooling function after completion of the program with opened firing chamber.

The quick cooling remains activated until the preheating temperature is reached.

In this case, the actual temperature of the firing chamber must be higher than the pre-heating temperature. Press the sensor button “☒” again to deactivate quick cooling.

Note: See section 6.2 for programmed quick cooling.

Press button “Start/Stop”

Press this button to start or stop a program. To start a program, the firing chamber must be open. Stopping the program is delayed by a few seconds. Please wait until the firing chamber is completely open and the furnace has returned to firing data view.

Clear “C”

Press this sensor button to:

- delete a wrong entry.
- confirm an informational or error display.
- abort a vacuum after reaching the preset vacuum level.
- abort a functional test.

4.5 Main menu

Figure 13: Main menu

Press the sensor button “M” to open the main menu. The main menu includes all functions that cannot be executed directly. The first line of the screen shows the menu name and the second line shows a sub-menu

selection or a parameter. Activate the individual sub-menus by pressing the sub-menu sensor buttons; activate presets in the same way. Change values with the numeric keypad. Changed values are saved by selecting a new sub-menu or pressing the sensor button “M” until the program is displayed again.

4.5.1 Configuration


Figure 14: Configuration menu

Sort by number/ name

Programs can either be sorted by their program number (numerically) or by the program name (alphabetically) (see section 6.12).

Note: This function can only be used with Latin characters.

Job number “on”

(inactive)

Used to enter job numbers in case of multiple firing runs.

Display time/date

Select whether time or date should be displayed in the program (see section 6.15).

Test mode

(only available in service mode)

Used to record the last software actions.

Furnace Press/Touch

(only available in service mode)

Used to select between Touch and Press functions.

Power failure backup duration 1–20 seconds

Used to adjust the backup time in case of power failure in seconds. Factory setting is 10 seconds (see section 6.4).

Date mode

Used to select between these date modes:

DD.MM.YYYY

YYYY.MM.DD

MM.DD.YYYY

Press the sensor field until the desired date mode is displayed.

Automatic activation on/off

Used to activate the automatic power-on time (see section 6.7).

Time 00:00 – 24:00

Used to set the time at which the furnace automatically switches on, if connected to main power. For safety reasons, this time can be set max. 24 hours in advance (see section 5.7).

4.5.2 Furnace parameters


Figure 15: Furnace parameters menu

Reset calibration offset, hours of operation and language settings here. All other displays in this menu are for informational purposes only and cannot be changed.

4.5.2.1 Reset calibration offset

Note: Reset the calibration to 0 °C prior to starting the silver test sequence.

1. Press sensor button “Calibration offset”.
2. Press sensor button “C”, numerical values are reset to zero.
3. Press sensor button “M” until the previous program is displayed.

Now open the silver test sequence (#475) and recalibrate the furnace (see section 6.5).

Note: The offset cannot be manually entered or changed. Therefore, always run the silver test sequence for calibration.

4.5.2.2 Vacuum test

1. Press sensor button “Vacuum test”.
2. Press sensor button “Start/Stop”.

The vacuum test starts.

Abort vacuum test by pressing the sensor button “Start/Stop”.

The set vacuum of 50 hPa must not decrease by more than 20 hPa within 5 minutes. After the program is complete, the screen displays whether this requirement is met.

4.5.2.3 Reset language

1. Press sensor button “Reset language”.
2. Confirm with the “C” button.
3. Disconnect **Multimat2 Touch** from main power.

4. Reconnect **Multimat2 Touch** to main power. The start screen is displayed, and the language must be selected. All other individual settings remain set.
5. . Select new language.

4.5.2.4 Heating % (power factor)

This function specifies, with which percentage of the main power the heating muffle operates. Depending on main power and heating muffle, the furnace automatically sets to 30 %, 75 %, 85 % or 100 %.

When ordering a replacement heating muffle, always specify the voltage on the type label of the furnace.

100; 110; 115; 125 V = 100 V muffle

230 V = 230 V muffle

4.5.2.5 Reset hours of operation (muffle)

This function can be used to reset the hours of operation of the muffle after replacing to zero.

1. Press the sensor button “Muffle hours”.
2. Enter the code “6070” with the numeric keypad.
3. Select “Muffle hours” again.
4. Confirm resetting with “C”.

4.5.2.6 Statistical data

Hours of operation (furnace)

Hours of operation (pump)

Lift cycles

Firing cycles

These values are read-only and are used as reference for maintenance cycles.

4.5.3 Smart Media Card

The Smart *Media* Card is used to save individual programs and for transfer to other **Multimat2 Touch** and **Multimat2 Touch+Press** furnaces. (Updating the operating system via Smart *Media* Card will soon be available, see also special functions, chapter 6.16.)

4.5.4 Delete programs

Figure 16: Program deletion menu

4.5.4.1

Single program

1. Press sensor button “M”.
2. Press sensor button “Delete programs”.
3. Press sensor button “Single program”.

4. Enter program number with the numeric keypad.
 5. Press sensor button “Single program” again.
- The single program is deleted.

4.5.4.2 All individual programs

1. Press sensor button “M”.
 2. Press sensor button “Delete programs”.
 3. Press sensor button “All individual programs”.
 4. Confirm the deletion by pressing “C”.
- All individual programs are deleted.

4.5.4.3 Delete fixed and special programs

The sensor fields of the fixed and special programs cannot be deleted.

5 Create individual programs

The individual steps until the ceramic selection have already been described in section 3.4 **First operation**.

5.1 Firing method


Figure 17: Firing modes menu

Dentsply DeTrey, the furnace changes to this screen:

After selecting a ceramic type, the firing methods are displayed. Select from this screen the desired firing method by pressing the appropriate sensor button. After selecting the desired firing method, e. g.

5.2 Program type

Figure 18: Program types menu

Select from the displayed possibilities by pressing the appropriate sensor button, e.g. individual programs. After selecting the desired program type, the furnace changes to the next screen:

5.3 Program list

Figure 19: Metal ceramics menu

Use the buttons “↑” and “↓” on the sensor bar on the bottom of the screen to scroll the program list. Press the button “NEW” to create a new program. The furnace changes to this screen:

5.4 Text input

Figure 20: Entering characters

Complete the entry by pressing “→”. The furnace changes to the next screen. Different characters are available, depending on the selected language.

5.5 Firing parameters


Figure 21: Form for entering firing parameters

Numerical values are entered with the numeric keypad.

Note:

For time data, the colon (:) between minutes:seconds and hours:minutes is a default setting and needs not to be

entered separately. Enter a time of 3 minutes and 20 seconds as 0320.

Note: The program name can be changed afterwards with the sensor button “Name”.

5.5.1 Set preheating temperature

Press the relevant yellow sensor button. Enter the temperature value with the numeric keypad, e.g. for 500 °C enter the value 500.

5.5.2 Set predrying time

(firing chamber in upper end position)

Only used for low-temperature drying for an extended period. During this drying phase, the firing chamber remains in the upper end position.

5.5.3 Set drying time

Press the relevant yellow sensor button. Enter the time value with the numeric keypad, e.g. for 6 minutes enter 0600 (see section 6.13).

5.5.4 Set preheating time

Press the relevant yellow sensor button. Enter the time value with the numeric keypad, e.g. for 3 minutes enter 0300.

5.5.5 Set vacuum level

Press the relevant yellow sensor button. Enter the vacuum level with the numeric keypad, e.g. for 50 hPa enter the value 50.

5.5.6 Set heating rate

Press the relevant yellow sensor button. Enter the heating rate with the numeric keypad, e.g. for 80 °C/min enter 80. In other firing modes, the heating rate can be time-controlled.

5.5.7 Set firing temperature

Press the relevant yellow sensor button. Enter the firing temperature with the numeric keypad, e.g. for 940 °C enter 940.

5.5.8 Set vacuum

a. time-controlled (e.g. Dentsply DeTrey mode)

Press the relevant yellow sensor button. Enter the vacuum time with the numeric keypad, e.g. for 1 minute enter 0100.

b. temperature-controlled (e.g. Dentsply Ceramco mode)

Press the relevant yellow sensor button. Enter the temperature for switching the vacuum on or off with the numeric keypad.

Note: The ceramic types “metal ceramics” and “pressable ceramics” are entered and displayed in minutes:seconds and the “sinter ceramics” in hours:minutes.

5.5.9 Set firing time

Press the relevant yellow sensor button. Enter the firing time with the numeric keypad, e.g. for 2 minutes enter 0200.

Note: The ceramic types “metal ceramics” and “pressable ceramics” are entered and displayed in minutes:seconds and the “sinter ceramics” in hours:minutes.

5.5.10 Set tempering temperature

During tempering, the thermal expansion coefficient (TEC) of the ceramic is increased by controlled growth of leucite crystals. Thus, the TEC of the metal ceramic can be adjusted to alloys that differ greatly in terms of their TEC.

Press the relevant yellow sensor button. Enter the tempering temperature with the numeric keypad, e.g. 1000 for 1000 °C.

Consult the ceramics manufacturer before using the tempering function!

5.5.11 Set tempering time

Press the relevant yellow sensor button. Enter the tempering time with the numeric keypad, e.g. enter 1000 for 10 minutes.

5.5.12 Set cooling stage

Using the 3 cooling stages of the **Multimat2 Touch** results in a stress relaxation in the ceramics. When cooling stages were programmed, the cooling starts after the end of the firing time, and lasts until the preheating temperature is reached again.

For activating it, press the yellow sensor button until the desired cooling stage is displayed.

Cooling stages

0 = Firing chamber moves immediately to upper end position – no controlled cooling

1 = Firing chamber opens about 70 mm

2 = Firing chamber opens about 50 mm

3 = Firing chamber remains closed

Quick cooling

(Controlled cooling is only available in sinter programs.)

5.6 Save program

Save your program by pressing the button “S”. The entered firing parameters will be checked on validity. Invalid firing parameters will be shown in the info field and must be corrected after pressing the button “C”. Once all firing parameters are accepted, the control system proposes to use the next available program number as a storage location. Press the button “S” again to accept the default setting.

If you wish to save a program under an individual number, use the numeric keypad to enter a different program number and confirm with “S”. The firing parameters will be saved. By entering a new number, the old program with this number will be overwritten. You will be asked to confirm the overwriting.

5.7 Start program

1. Open the firing chamber with the button “↑” on the right sensor keypad (if not opened yet).
2. Press the sensor button “Start/Stop”. The furnace changes to the screen “Firing curve”.

The firing curve illustrates all firing steps. Initially, the firing curve is empty, and it will be filled with color during the firing sequence. Below the firing curve the times for the individual firing steps are displayed, counting down to zero. The displayed remaining time is the estimated addition of all partial times.

If selected in the basic settings, the printer or PC icon is displayed on the right edge of the screen.

The following functions can be selected during firing:

Function	Activation
a. View	Sensor button lift „↑“, and „↓“ to close again
b. Abort view	Sensor button “Start/Stop”
c. Abort program	Sensor button “Start/Stop”
d. Abort vacuum/view	Sensor button “C”
e. Change firing parameters	Sensor button “←”
f. Confirm messages	Sensor button “C”
g. Activate Night Mode	Sensor button “Night Mode”

Programming the firing parameters is the same in all firing modes.

5.8 Standard functions

5.8.1 Opening a program

1. Press the sensor button “R” and
2. enter the program number.
3. For program numbers lower than 100, press “R” again.
4. Press the soft key “List” and browse the available programs. The desired program is directly displayed by pressing the relevant sensor button.

5.8.2 Saving a program

Save your program by pressing the button “S”. The entered firing parameters will be checked on validity. Invalid firing parameters will be shown in the info field and must be corrected after pressing the button “C”. Once all firing parameters are accepted, the control system proposes to use the next available program number as a storage location. Press the button “S” again to accept the default setting.

If you wish to save a program under an individual number, use the numeric keypad to enter a different program number and confirm with “S”. The firing parameters will be saved. By entering a new number, the old program with this number will be overwritten. You will be asked to confirm the overwriting.

5.8.3 Changing a program

Change your program by directly selecting the firing parameter and entering the new value with the numeric keypad.

The new value will be saved by:

- a. selecting the firing parameter again or
- b. selecting another firing parameter or
- c. saving the program

5.8.4 Changing a program during the firing sequence

All firing sections not yet processed may be changed while firing is in progress. For this, press the button “←” in the firing curve screen. The firing parameter screen is displayed. The firing process can also be observed in this screen. Select the firing parameter to be changed, enter changed value and confirm by selecting the parameter again (see section 5.8.3). The yellow sensor buttons of the processed and running firing section are masked out and cannot be changed. The change is one-time, is not saved and is only valid for the currently running program. Press the button “⇒” to change to the firing curve screen again. After the firing is complete, the original parameters are displayed again.

5.8.5 Copying a program

Open the program to be copied, Press the sensor button “S” and enter a new program number with the numeric keypad, if required. Press the button “S” again to copy the program. If the entered program number already exists, a message is displayed.

Press the sensor button “S” twice to copy the program to the next free program number.

5.8.6 Deleting programs

See section 4.5.4

5.9 Pressing

In order to press all-ceramics, first place the compression platform on the firing socket holder.

After selecting the ceramic type "Pressable ceramic" in the screen "Ceramic type", the furnace changes directly to the program list "Pressable ceramics". The compression program for FINESSE ALL CERAMIC is saved as fixed program at the program number 301.

Select the factory-preset program FINESSE ALL CERAMIC from the program list by pressing the appropriate sensor button. After selecting the desired program, the furnace automatically changes to the screen FINESSE ALL CERAMIC.

Open the firing chamber with the sensor button "↑" and start the pressing program with the sensor button "Start/Stop". The furnace changes to the Program procedure screen.

Figure 22: Pressing program menu

Note: Do not place the investment ring into the pressing muffle at this time.

Press procedure description:

- After the start, the firing chamber closes and the temperature increases from idle temperature (adjustable in the basic settings) to the starting temperature of 700 °C.
- The firing chamber opens automatically as soon as the starting temperature has been reached. The upper end position of the firing chamber is 158 mm, as to facilitate placing the pressing muffle on the pressing socket. An acoustic signal informs the user to place the externally preheated pressing muffle onto the pressing socket
- After placing the pressing muffle, press the sensor button Lift "↓", and the pressing program continues.
- The firing chamber closes; the vacuum pump activates and generates a vacuum of 50 hPa.
- After reaching the vacuum level, the temperature increases with a heating rate of 60 °C/min to the pressing temperature of 930 °C.
- The hold time at pressing temperature is 20 min.

- After completion of the hold time, the pressing cylinder is loaded with 2.7 bar compressed air and the pressing procedure starts. The compression force is 21.7 dN, and the pressing time 10 minutes.
- At the end of the pressing time, the pressing cylinder moves into its initial position
- The pressing chamber is ventilated and moves to its upper end position. The pressing is now complete.

Note: Always observe the manufacturer's requirements when processing pressable ceramics. Always make sure that the pressing muffle is centered on the pressing socket.

Caution! Use the pressing socket only for pressing. Crowns and bridges fired on the pressing socket will be underfired.

6 Special functions

6.1 Night/Standby

This function prevents moisture from entering the firing chamber by maintaining a preselected temperature between 101 °C and 300 °C in the firing chamber when switched "Off". The Night/Standby temperature can be selected in the basic settings.

Recommended Night/Standby temperature: 120 °C.

Temperature > 100 °C = Night/Standby activated

Temperature < 100 °C = Night/Standby deactivated

Activating the Night/Standby function

- Night/Standby temperature greater than 100 °C and switching the furnace off with the "On/Off" switch.
- Night/Standby temperature greater than 100 °C and pressing the sensor button with the bed symbol during firing (Night mode).

Note: In order for this function to be active, the furnace must not be switched off (e.g. main switch/master switch must not be switched off).

6.2 Quick cooling

The function quick cooling is used to decrease the temperature of the firing chamber faster after firing is complete. On activation, the firing chamber opens, the vacuum pump activates and continuously aspirates fresh air into the firing chamber. Quick cooling ends when the preheating temperature has been reached. This means that the pump will stop 20 °C below the preheating temperature.

Manual activation:

Quick cooling can be manually activated with opened firing chamber after completion of the program by pressing the sensor button "■".

In this case, the actual temperature of the firing chamber must be greater than the preheating temperature.

Switch quick cooling off by pressing the sensor button “” again.


Automatic activation:

If quick cooling should be activated automatically with the selected program, you have to press the sensor button “Cooling stage” in the firing parameter screen until quick cooling is displayed.

The symbol for quick cooling “” is displayed in the firing curve screen.

6.3 View function

(for soldering)

During soldering without vacuum it is permissible to move the firing chamber upwards with the lift button “” and stop it with the lift button “”, as to inspect the soldering material. The soldering time (firing time) will be interrupted, and the screen now switches automatically to the firing data screen. You can now change the soldering (firing) temperature.

1. Press sensor button “Firing temperature”.
2. Enter the new firing temperature code with the numeric keypad
3. Press sensor button “Firing temperature” again.
4. Press the sensor button “”. The screen changes to the program sequence screen.
5. Press the sensor button Lift “”.

The firing chamber closes again and the program continues. After the temperature has been set, post-firing is possible for up to three minutes. The post-firing time counts up from zero. Post firing can be aborted at any time with the “Start/Stop” button.

Note: The View function can only be activated when tempering or cooling is not activated. The temperature change is only visible in the actual value. The theoretical value remains unchanged.

6.4 Power failure backup

(sub-menu configuration section 4.5.1)

The **Multimat2 Touch** is equipped with an adjustable power loss backup. This backup activates as soon as there is a main power failure during a running firing program. The backup duration can be set from 1 – 20 seconds. If the downtime is shorter than the backup time, the program will continue to run. At the end of the program the following message appears: “E 07 – There was a power failure”.

If the downtime is longer than the backup time, the program will be aborted and the following message is also displayed: “E 07 – There was a power failure”.

Note: Always check the firing result!

Adjust power failure backup

1. Press sensor button “M”.
2. Press sensor button “Configuration”.
3. Press sensor button “Backup time”.
4. Set backup time with the numeric keypad between 1 and 20 seconds.
5. Press sensor button “M” until the previous program is displayed.

6.5 Testing the firing chamber temperature

(with silver calibration set “silver sample manual”, REF D03 532 803)

The firing chamber temperature is precisely set in the factory. Should checking the firing chamber temperature become necessary for any reason, please proceed as follows:

Required items:

- 1 Dentsply silver sample carrier
- 1 piece of silver wire – diameter: 0.3 mm, length: 37 mm

Procedure:

- Allow furnace to heat through for approx. 1 hour at 600 °C.
- Place silver wire into the Dentsply silver sample carrier.
- Reset calibration offset (see section 4.5.2.1).

The following data are saved as special programs under “silver sample”, program number 375:

Preheating temperature	650 °C
Predrying	0 minutes
Drying	0 minutes
Preheating	3 minutes
Vacuum level	0 hPa
Heating rate	120 °C/min
Firing temperature	961 °C (melting point of the silver wire)
Vacuum time	0 minutes
Firing time	1 minute
Tempering temperature	0 °C
Tempering	0 minutes
Cooling stage	0

1. Open firing chamber and place sample carrier with silver wire centrally onto the firing platform.
2. Start test program “Silver wire calibration”.
3. If after completion of the program the silver wire has started to melt along its upper third, i.e. a bead has formed on the wire surface, then the temperature is correct (with an accuracy of ± 2 °C).
4. If the silver wire surface has not started to melt, then the firing chamber temperature is too low. In this case, repeat the test program – each time with a temperature increase of 3 °C – until the required melting effect on the surface of the silver wire is achieved.
5. If the silver wire has melted altogether, then the firing chamber temperature is too high. In this case, repeat the test program – each time with a new length of silver wire and a temperature reduction of 3 °C – until the required melt effect on the surface of the silver wire is achieved.

6. After you have found the correct melting temperature, press the sensor button “Silver”. The temperature correction is taken over as “Calibration offset”.

The correction of the firing chamber temperature is completed.

6.6 Night Mode

This function switches the heating off automatically after firing is complete and the chamber temperature drops to the night/standby temperature. This function can be switched on or off for the program as desired by pressing the “Bed” sensor button. The “Bed” symbol for night mode is displayed in the firing curve. After firing is complete, the furnace switches the display off and the firing chamber cools down. After reaching the night/standby temperature, the firing chamber closes automatically. The night/standby temperature is maintained until the furnace is switched on again. Pressing the “On/Off” sensor button switches the **Multimat2 Touch** on again.

6.7 Select activation time

(sub-menu Configuration section 4.5.1)

This function automatically activates the furnace at a preselected time. The function is set in the sub-menu “Configuration”. The input fields “Activate” and “Time” are displayed.

Setting the activation time

1. Press sensor button “M”.
2. Press sensor button “Configuration”.
3. Press sensor button “Time”.
4. Enter the 24-hour time in the form 00:00 with the numeric keypad.
5. Set sensor button “Activate” to “on”.
6. Press sensor button “M” until the previous program is displayed.

When switching the furnace off with the “On/Off” switch, it will be automatically activated again at the preset time.

Note: For safety reasons, this function is active for only 24 hours. Therefore the sensor button “Activate” will switch back to “Off” after automatic activation and must be manually set to “On” again.

6.8 Vacuum test

(sub-menu Furnace parameters, section 4.5.2.2)

This test is for service purposes only.

6.9 Print/PC

Printing

The theoretical and actual values of a firing program are automatically sent to a printer or PC after completion of the firing, if the data output is set to “Printer” or “PC” in the basic settings. The theoretical values for firing data can be printed out at any time by pressing the button PR/PC. If no printer or PC symbol is displayed, the data output is disabled.

The printer must be configured as follows (please see your printer manual for directions):

9600 baud / 8 bit / 1 stop bit / no parity

Note: For printer/PC connection you will require:

- for serial connection:
1 “twin” cable, serial (D-SUB, 9, St/D-SUB, 9 jacks, 9 wires)
or
- for parallel connection:
1 serial to parallel connector, if a parallel printer is available.

Connect the printer with the **Multimat2 Touch** with the RS 232 interface.

6.10 Acoustic signals

Short signal:	At each key press
Long signal:	At invalid inputs
Triple signal:	At program end
	In the pressing program: insert material/furnace is preheated

6.11 Software update

In the scope of ongoing product maintenance, we recommend always installing the newest software version for the **Multimat2 Touch** and **Multimat2 Touch+Press** products. Please ask your field agent or simply download the newest software version from the internet. The internet address is: www.digitux.de/dentsply. Password: “tech”.

The software of the **Multimat2 Touch** or **Multimat2 Touch+Press** is updated via the serial interface. Therefore, please connect your Multimat2 with a serial cable to your PC.

1. Boot PC / Laptop
2. Open the Windows Explorer (Start/Programs/Windows Explorer).
3. Insert floppy disc or CD-ROM containing the new software.
4. Click on 3.5" Disc (A:) or CD-ROM.
5. On the right hand side of this window you should see a file called MMT_Updt.exe, which has a satellite dish as a symbol.
6. Click on this file.
7. Choose your language from the options shown.
8. Disconnect the furnace from the mains supply (pull out the plug).
9. Re-connect the furnace to the mains supply (put the plug back in the socket).
10. The computer will now start with the data transfer (approx. 2–3 minutes). The update is finished when "Ready!" appears in the dialogue window of your furnace.
11. Disconnect the power cable and the serial cable from the furnace. Wait for 5 seconds, and then reconnect only the power cable.
12. The **Multimat2 Touch** display now shows the latest software version on the start screen.

6.12 Program sorting by name or program number

(Main menu, section 4.5.1)

This function is used to sort programs by name or program number. As a helpful feature, programs are given a letter prefix for sorting alphabetically.

Program sorting

1. Press sensor button "M".
2. Press sensor button "Configuration".
3. Press sensor button "Sorting" and select by name or number.
4. Press sensor button "M" until the previous program is displayed.

Note: This function can only be used with Latin characters.

6.13 Change standard drying position and steps

This function can be used to change the standard drying position and the number of steps after which the standard drying position is reached, measured from the upper end position.

Set new drying position

1. Press the text field "Drying" (not the sensor button!) twice. A selection screen is displayed.
2. Press the sensor button "Drying position".
3. Enter the new drying position in "mm" with the numeric sensor keypad.
4. If the drying position is lower than 100 mm, press the sensor button "Drying position" again.

Note: The drying position can be adjusted between 0 and 150 mm. 0 mm is the starting position in the selection screen.

Adjust drying steps

1. Press the sensor button "Steps".
2. Enter number of steps with the numeric keypad.
3. Press "Arrow left", and the program is displayed again.

Note: To permanently store the changed drying parameters, the program must be saved before the start. The new drying parameters apply only for the changed program. All other programs remain unaffected.

Caution! We recommend not changing the default standard drying positions and steps without urgent reasons.

6.14 Change tempering positions

This function can be used to change the tempering position.

Set new tempering position

1. Press the text field "Tempering" in the program twice. A selection screen is displayed.
2. Press sensor button "Tempering position".
3. Enter the new tempering position in "mm" with the numeric sensor keypad.
4. If the tempering position is lower than 100 mm, press the sensor button "Tempering position" again.
5. Press "Arrow left", and the program is displayed again.

Note: To permanently store the changed tempering position, the program must be saved before the start. The new tempering position applies only for the changed program. All other programs remain unaffected.

**Change the tempering positions for
ceramics only after consultation with the
manufacturer.**

6.15 Display date or time

This function is used to set either the date or the time to be displayed in the program. The factory setting is the date.

Changing the time or date display

1. Press the sensor button "Drying position".
2. Press the sensor button "Configuration".
3. Press the sensor button "Display date" or "Display time". The display changes to date or time.
4. Touch sensor button "M" until the previous program is displayed.

6.16 Smart Media Card

(Smart Media Card REF 54 6023 0100)

(Main menu section 4.5)

Data are backed up on a Smart Media Card (see Figures 23 and 24). This card is a dedicated storage medium, available at capacities of up to 128 MB. The cards come in two types:

- Cards for an operating voltage of +5 V, characterized by its left upper corner being cut off.
- Cards for an operating voltage of +3.3 V, characterized by its right upper corner being cut off (see Figure 23). In addition, this card type is usually labelled with the letters "ID".

Either card at any capacity can be used for the Multimat 2 unit. Both types of cards feature a write-protection field. If this field is covered by an adhesive conductive metal foil, the Smart Media Card will be write-protected (see Figure 24).


Figure 23: +3.3 V Smart Media Card with 128 MB of memory


Figure 24: +3.3 V Smart Media Card with write protection in place

Insert the Smart Media Card into the card reader with the medal surface facing back (outward, toward you) (see Figure 25).

To remove the card, press the eject button (see Figure 26).

Using the Smart Media Card function, all or individual functions of the **Multimat2 Touch** can be saved to the Smart Media Card and restored from the Smart Media Card to the **Multimat2 Touch**. The Smart Media Card, a dynamic storage medium, can be used to save several hundred individual programs. This is particularly useful for additional program backup.

Figure 25: Inserting the Smart Media Card


Figure 26: Smart Media Card eject button

Format Smart Media Card

Figure 27: Smart Media Card menu

Prior to first use, the Smart Media Card must be formatted.

In doing so, all data on the Smart Media Card will be deleted.

1. Insert the Smart Media Card into the slot on the right rear side of the furnace (metal surface facing back).
2. Press the sensor button "M" on the numeric keypad.
3. Select "Smart Card".
4. Select "Format Smart Card".
5. Confirm with "S" = yes (abort with "C" = no).
6. Formatting will take about 5 to 10 seconds.
7. Complete successful formatting by pressing the "C" button.
8. Now the individual programs can be stored on the card.

6.16.1 Program transfer furnace → Smart Media Card

1. Insert the Smart Media Card into the slot on the right rear side of the furnace (metal surface facing back).
2. Press the sensor button “M” on the numeric keypad.
3. Select “Smart Card”.
4. Press the sensor button “Furnace → Smart Card”.
5. Select individual programs from the list or press “All” in the lower soft key bar.
6. Press “→ Card” in the lower soft key bar.
7. Confirm program saving by pressing “C”.

6.16.2 Program transfer Smart Media Card → furnace

Using this function, you can transfer programs from the Smart Media Card to the **Multimat2 Touch**.

1. Insert the Smart Media Card into the slot on the right rear side of the furnace (metal surface facing back).
2. Press the sensor button “M” on the numeric keypad.
3. Select “Smart Card”.
4. Press the sensor button “Smart Card → Furnace”.
5. Confirm loading programs with “S” = yes (abort with “C” = no).
6. The program(s) is/are being loaded.
7. Confirm program saving with the “C” button.

Service note:

In the event that your firing programs for older Multimat models were previously stored on a Smart Card (see Figure 28) and you want to continue using your data on a Multimat2, you can have the data transferred to a Smart Media Card, free of charge, by your distributor.

Figure 28: Smart Card, front and rear view

7 Service and maintenance

7.1 Lift

The lift is used for vertical transport of the firing chamber. The unit is driven by a geared motor. If there is a power failure, the firing chamber can be lifted upwards by hand in order to remove the work pieces. Apply vigorous pressure downwards to move the firing chamber downwards by hand.

Caution! Do not change the setting of the belt tension screw. A too strong belt setting may cause overload of the lift drive, thus damaging the unit.

7.2 Firing socket/pressing socket and Firing socket carrier

The firing platform carrier is mounted vertically resilient on the top of the furnace subsection. The O-ring in the firing socket carrier seals the firing chamber during vacuum firing. Depending on whether the unit fires or presses, the firing socket or the pressing socket is placed on the firing socket carrier. Both firing and pressing socket insulate the firing chamber against downward temperature loss and are also used as object carrier.

Maintenance note: The O-ring of the firing socket carrier must be kept clean and should be occasionally checked for damage.

Caution! In order to prevent damage to the firing socket carrier and the O-ring seal, either the firing socket or the pressing socket must be placed on the firing socket carrier.

7.3 Vacuum pump

We recommend operating the **Multimat 2 Touch** with the highperformance Dentsply vacuum pump. However, the **Multimat2 Touch** can be operated with any highperformance vacuum pump with grounded plug connector (power uptake max. 2.5 A) that has a liter performance of at least 30 l/min and a final pressure of 30 hPa. Connecting the pump is described in section 3 “Setup and first operation”. The pump connection tube must not be longer than 2 meters

Maintenance note: Observe the operating manual of the pump for maintenance purposes.

Important: Vacuum pumps with oil lubrication must have an oil change every 3 months.

7.4 Changing the firing muffle

(trained electricians or service personnel only!)

Caution!

This product contains ceramic fibers, and can release fiber dusts. Fiber dusts have proven to be carcinogenic in animal experiments. Observe EC safety data sheet.

The heat insulation of the firing chamber in the **Multimat2 Touch** consists of ceramic fibers. After extended use of ceramic fibers at temperatures above 900 °C, silicone type substances (cristobalite) may be generated. In certain cases, e. g. when replacing the firing muffle, dust exposure may occur, which may cause irritation of the skin, eyes, and respiratory organs. Proceed as follows when replacing the firing muffle:

- All personnel must wear long sleeved clothing, headwear, eye protection and gloves.
- Attach dust exhaustion to the dust source; if not possible, equip personnel with dust mask FFP3 or similar.
- After completion of the tasks, possible adhering dust must be first washed off the unprotected skin. Afterwards wash with warm water and soap.
- Wash working clothes separately.

Warning!

Disconnect the Multimat Touch from main power before opening the housing.

Service work may only be carried out by qualified professionals or service technicians.

Disassembly:

Multimat2 Touch+Press

1. *Observe dust protection measures!*
2. *Disconnect from main power!*
3. *Switch off compressed air feed!*
4. *Unscrew pressing hood, pull backwards somewhat and carefully remove upwards.*
5. *Disconnect plug connector from the solenoid valve.*
6. *Release the compressed air tube from the solenoid valve by pushing the clamping ring upwards and simultaneously pulling off the compressed air tube.*
7. *Disconnect final position sensor.*
8. *Proceed further as described below for the **Multimat2 Touch**.*

Multimat2 Touch

1. Unscrew cover and lift off vertically upwards (until pressing punch is free).
2. Disconnect protective conductor from the cover.
3. Remove insulating disc(s).
4. Disconnect thermal element from both rear connecting bolts.
5. First remove the thermal element together with the bar, and then the terminal insulation.
6. Disconnect the heating wire ends of the muffle from the two front connecting bars (use a 12 mm flat wrench for the anti-twist device of the connecting bar and a 7 mm flat wrench for loosening the hexagon bolt).
7. Remove muffle from the insulating insert; be careful if the furnace is hot and touch only the ends of the heating wire.

Assembly:

Caution!

Please wear cotton gloves when installing the new muffle, as to prevent sweat from the hands on the quartz tube.

Installation of a new muffle is done in reverse order of the disassembly.

Note: After a muffle change, we recommend testing the temperature with the silver test (see section 6.5).

Warning!

Check that the wire ends are not in contact with the metal mantle of the upper furnace section.

Reconnect protective conductor to the cover! Do not shut the protective conductor between vessel and cover.

7.5 Replacing the control unit

(trained electricians or service personnel only!)

Warning!

Disconnect the unit from main power before opening the housing.

Service work may only be carried out by qualified professionals or service technicians.

Caution!

The firing chamber must be in its upper end position.

1. Disconnect from main power.
2. Unscrew all Phillips screws below the control element.
3. Remove control unit in upward direction.
4. Disconnect vacuum tube and all plug connections.
5. Disconnect thermal sensor.
6. Disconnect protective conductor from the control unit.

Installation of the new control unit is done in reverse order.

Note: For shipping the control, use the reusable packaging only!

7.6 Inspections

After 15,000 lift cycles a warning note is indicated in the display which can be reset with the "C" button. We recommend an inspection of the furnace. During this inspection, essential functions are checked and components replaced, if required. Please ask your specialist dealer for Dentsply services.

7.7 Cleaning notes

The **Multimat2 Touch** is varnished with an epoxy paint. Please use a mild household detergent for cleaning the varnished surface. The touch screen may be cleaned only when the unit is switched off with a mild window cleaner.

8 Faults and remedies

The following information is intended to assist you in identifying and evaluating simple errors and to fix them by appropriate remedial action.

An error message, its probable cause and the appropriate remedial action is shown in brief on the display. Confirm the displayed message by pressing the “C” key. Additional information to the reason of the fault can be displayed by pressing the “?” button in the left upper corner.

For faster identification, the messages are labeled with an alpha-numeric prefix. If the fault cannot be corrected, please contact your distributor. Spare parts can be ordered by your distributor only.

Warning!

Disconnect the unit from main power before opening the housing.

Caution!

It is the users' responsibility to send the furnace back in a proper package.

Use only the original package. Contact the supplier if you need a new original package.

Confirm all error and information messages first with “C”.

No.	Message	Cause	Remedy
	Screen is black and green power LED does not light up.	Display faulty. Power fuses faulty.	Disconnect from main power! Replace faulty fuses. If the display remains black, contact service technician.
E 01	Waiting time for inserting object exceeded.	The 15 minute waiting time for inserting the pressing muffle was exceeded.	Expedite insertion.
E 02	Vacuum is not removed, firing chamber does not open.	Solenoid valve faulty.	Ventilation valve faulty. Contact service technician.
E 03	Limit switch not closed.	Transport protection or foreign body stuck between upper and lower part of furnace.	Remove transport protection or foreign body.
E 04	Control unit is overheated. Program aborts, continuous acoustic signal.	Extremely high standby temperature with open firing chamber.	Remove power plug. Allow furnace to cool for approx. 5 minutes. Restart furnace after the cooling time.
E 05	Theoretical vacuum not reached.	Vacuum system leaks. Suction power of the pump is too low. Moisture inside the firing chamber insulation.	Start vacuum test program (in the menu unit parameters). Observe pressure display, close supply via the vacuum tube (bend tube) after 1 minute. If the vacuum drops sharply, the firing chamber leaks or there is moisture in the firing chamber. If the vacuum does not drop, the pump has insufficient suction power. Check pump according to manufacturer's instructions or use stronger pump. Start dehumidifying program 476 with vacuum pump or 474 without vacuum pump.
E 06	Quick cooling is not complete.	Preheating temperature is not reached yet.	Wait until idle temperature is reached or abort with "C".
E 07	There was a power failure.	Temporary power failure during firing.	Please check the firing result!
E 08	Heating circuit is faulty. Program aborts.	Muffle or Triac faulty.	Contact service technician.
E 09	Thermal element failure or thermal wiring failure. Program aborts.	a. Thermal element broken or thermal wiring disconnected. b. Thermal element poles reversed.	a. Contact service technician. b. Reverse polarity of thermal element or thermal wiring.
E 10	Battery weak.	Battery on the control board empty.	Have service technician replace battery.
E 11	Overheating.	The actual temperature is 25 °C above the theoretical value.	Have service technician check the temperature control.
E 12	Maximum pressing stroke reached. Program terminated.	Insufficient ceramic material in the muffle. Sensor not positioned or faulty.	Use sufficient amount of ceramic material. Contact service technician.
E 13	Sensor for end position of pressing not detected.	End position sensor misaligned or faulty	Contact service technician.
E 17	Temperature reference generator faulty.	Temperature reference generator faulty!	Contact service technician.
E 20	Smart <i>Media</i> Card cannot be read.	Smart <i>Media</i> Card inserted incorrectly. Reader faulty.	Insert Smart <i>Media</i> Card correctly. Contact service technician.
E 21	Smart <i>Media</i> Card is empty.	Empty Smart Card without data inserted!	Replace Smart <i>Media</i> Card.
E 22	Smart <i>Media</i> Card contains invalid data.	Smart <i>Media</i> Card was incorrectly written.	Write Smart <i>Media</i> Card again.
E 23	Not enough memory capacity in the furnace.	Not enough memory free in the furnace to save further programs!	Delete unused programs.

9 Preset programs

9.1 Firing table for Ceramco 3 metal ceramics (Dentsply mode)

(Use honeycomb firing material carrier.)

Prog. No.	Firing program	Preheat temp. (°C)	Drying (min)	Preheat (min)	Vacuum level (hPa)	Heating rate (°C/min)	Firing temp. (°C)	Vacuum time (min)	Firing time (min)	Cool stage
302	D-C3- Paste opaque	500	05:00	03:00	50	100	975	00:01	00:01	0
303	D-C3- Powder opaque	650	05:00	05:00	50	70	965	00:01	00:01	0
304	D-C3- Margin	650	05:00	05:00	50	70	965	00:01	00:01	0
305	D-C3- Opaque dentine/ Dentine effect/ Dentine enamel	650	05:00	05:00	50	45	930	00:01	01:00	0
306	D-C3-Glaze without Glaze liquid	650	03:00	03:00	0	45	920	00:00	03:00	0
307	D-C3- Glaze with Glaze liquid	650	03:00	03:00	0	55	925	00:00	00:30	0
308	D-C3- Correction (Add-On)	650	05:00	05:00	50	55	920	00:01	00:01	0

D = Dentsply mode

C3 = Dentsply Ceramco 3-metal ceramic

Firing temperatures are recommendations.

Make a test firing and adjust firing temperatures and times, if required.

For a large number of objects, temperatures or times must be increased appropriately.

9.2 Firing table for Finesse metal ceramic (Dentsply Mode)

(Use honeycomb firing material carrier.)

Prog #	Program Description	Preheat/ Low Temp. °C	Drying Time min.	Preheat min.	Vacuum Level hPa	Rate of Temp. Rise °C/min.	Firing/ High Temp. °C	Vacuum Time min.	Firing Time min.
311	D-FIN-Powder- Opaque	450	03:00	03:00	50	90	800	00:30	01:00
312	D-FIN-Paste- Opaque	450	05:00	03:00	50	90	790	00:06	00:30
313	D-FIN-Joint Opaque/ Margin/margin mod.	675	03:00	07:00	50	35	770	00:06	00:30
314	D-FIN-Correction Opaque	450	05:00	03:00	50	90	760	00:06	00:30
315	D-FIN-1 st Main Firing Opaque dentine, Dentine, Dentine effect, enamel, gums	450	05:00	05:00	50	35	760	00:06	00:30
316	D-FIN-2 nd Main Firing Opaque dentine, Dentine, Dentine effect, enamel, gums	450	05:00	05:00	50	35	750	00:06	00:30
317	D-FIN-Correction Firing	450	05:00	05:00	50	55	730	00:06	00:30
318	D-FIN-Glaze firing/ Stains*	450	03:00	03:00	–	70	750	–	00:06
319	D-FAC-repair porcelain	450	05:00	05:00	50	35	810	00:06	00:30
320	D-FAC-Shade stain porcelain	450	03:00	01:00	–	100	730	–	01:00

D = Dentsply mode

FIN = Finesse Metal ceramic

FAC = Finesse All Ceramic

* depending on desired degree of glaze

Caution!

Firing parameters not listed in the firing table must be set to “0” when entering values.

Caution! Do not use pressing socket for firing.

Notes: Make a test firing if required, and adjust firing temperatures or times. For a large number of objects, temperatures or times must be increased appropriately.

9.3 Firing table for Finesse All Ceramic (Dentsply Mode)

Caution! Use the supplied pressing socket.

Prog. No.	Pressing	Start temp. °C	Vacuum level hPa	Heating rate °C/min	Pressing temp. °C	Vacuum time min	Holding time min	Pressing time min
301	D-FAC-Press	700	50	60	930	27:00	20:00	7:00
319	D-FAC-repair porcelain	See firing table for Finesse Metal Ceramic / Progr. No. 319						
320	D-FAC-Shade stain porcelain	See firing table for Finesse Metal Ceramic / Progr. No. 320						

D = Dentsply mode

FAC = Finesse All Ceramic

9.4 Firing table for Finesse, FAC, Ceramco II (Dentsply Ceramco Mode)

(Use honeycomb firing material carrier as well as for pressing the supplied pressing socket.)

Prog #	Program Description	Time (min)					Temperature					Vacuum	Stages
		Pre-dry	Dry (min)	Pre-heat (min)	Vacuum hold	Firing time	Preheat/ Low Temp. (°C)	Firing/ High Temp. (°C)	Heate rate/ Rate of temp. rise (°C/min)	Vacuum on	Vacuum off	Vacuum level hPa	Cooling stages
324	C-FIN-Powder Opaque	0	03:00	03:00	0	01:00	450	800	90	450	800	50	0
325	C-FIN-Paste Opaque	0	05:00	03:00	0	0	450	790	90	450	790	50	0
326	C-FIN-Opaque/Correction/Opac	0	05:00	05:00	0	00:30	450	760	35	450	760	50	0
327	C-FIN-1st Dentin/Enamel/Mod.	0	05:00	05:00	0	00:30	450	760	35	450	760	50	0
328	C-FIN-2nd Dentin/Enamel/Mod.	0	05:00	05:00	0	0	450	750	35	450	750	50	0
329	C-FIN-Dentin Correction	0	05:00	05:00	0	0	450	730	55	450	710	50	0
330	C-FIN-Margin/Margin/Mod.	0	03:00	07:00	0	0	675	770	35	675	770	50	0
331	C-FIN-Natural Glaze	0	03:00	03:00	0	0	450	750	70				0
332	C-FIN-Glaze/Stain	0	03:00	03:00	0	0	450	750	70				0
333	C-CII-Ultra-Pake	0	03:00	03:00	0	01:00	500	975	200	500	975	50	0
334	C-CII-Paint-O-Pake	0	03:00	03:00	0	0	650	970	70	650	950	50	0
335	C-CII-Edge Shoulder	0	05:00	05:00	0	0	650	965	70	650	945	50	0
336	C-CII-Body	0	05:00	05:00	0	0	650	940	70	650	920	50	0
337	C-CII-Natural Glaze	0	03:00	03:00	0	01:00	650	940	70				0
338	C-CII-Low Temp. Glaze	0	03:00	03:00	0	01:00	650	915	70				0
339	C-CII-Add-On	0	05:00	05:00	0	0	650	915	70	650	865	50	0
340	C-CII-Silver Body	0	05:00	05:00	0	0	650	960	55	650	940	50	0
341	C-CII-Silver Natural Glaze	0	03:00	03:00	0	01:00	650	960	55				0
342	C-CII-Silver Overglaze	0	03:00	03:00	0	01:00	650	935	70				0
343	C-CII-Silver Add-On	0	05:00	05:00	0	0	650	940	70	650	890	50	0
344	C-CII-Colorlogic Veneer	0	08:00	08:00	0	0	500	940	50	500	920	50	1
345	C-CII-Colorlogic RDM – single	0	06:00	06:00	0	10:00	650	1010	50				1
346	C-CII-Colorlogic RDM – model	0	06:00	06:00	0	15:00	650	1010	50				2
347	C-CII-Final Touch Body & Incisal	0	05:00	05:00	0	0	450	705	55	450	685	50	0
348	C-CII-Final Touch Stains & Glaze	0	03:00	03:00	0	01:30	450	690	55				0
349	C-FAC-Repair Porcelain	0	05:00	05:00	00:06	00:30	450	770	35	450	770	50	0
350	C-FAC-stains	0	03:00	01:00	0	01:00	450	730	100	0	0	0	0

C = Dentsply Ceramco mode

FIN = Finesse metal ceramic

FAC = Finesse All Ceramic

CII = Ceramco II metal ceramic

10 Guarantee

Dentsply Furnaces

Section 1 Extent of the guarantee

1. The Manufacturer/Guarantor grants the Buyer/Beneficiary a Guarantee encompassing the functionality of the components listed in Section 2 for a period of two (2) years.
2. If one of the aforementioned components becomes defective during the warranty period, the Buyer/Beneficiary can demand to have the function restored professionally by replacement or repair of the defective component. This guarantee does not give the beneficiary the right to rescind the purchase agreement or to reduce the purchase price. The Manufacturer shall repair, at no charge to the Beneficiary, any defect of the Device covered by this Guarantee caused by a defect in material or workmanship, provided such defects are reported immediately after discovery.
3. The repair work covered by this Guarantee shall include any testing, measurements and adjustments required to repair a defect covered by this Guarantee; it shall not cover standard maintenance tasks recommended by the Manufacturer.

Section 2 Extent of the guarantee

1. The Guarantees shall include all component groups and components of the Device.
2. To be eligible for any claims under this Guarantee, the Beneficiary shall report the defect, at the same time presenting a copy of the original invoice.
3. Services under this Guarantee shall not extend nor re-start the term of this Guarantee.
4. The guarantee term for any replacement parts installed in the Device shall end concurrently with the guarantee term of the entire device.

Section 3 Exclusions

1. This Guarantee shall be void if any persons not expressly authorized by the Manufacturer/Guarantor make repairs or other modifications that are beyond the normal use of the Device or if a Device made by the Manufacturer/Guarantor is equipped with replacement parts not specifically approved for that Device.
2. This Guarantee shall not cover any transport damage, i.e. any damage caused by direct physical external mechanical force.
3. This Warranty shall be void if the Device is used improperly, in a way not described in the Instructions for Use, or for any other undue purposes.

4. The Manufacturer/Guarantor shall not be liable for any claims beyond those described in Section 1, including but not limited to liability for any damage caused outside the Device, unless liability is irrevocably assigned to the Manufacturer/Guarantor under the law.

Section 4 Obligations of the Buyer/Beneficiary

1. The Buyer shall have all maintenance, calibration and care tasks recommended for the Device in its Instructions for Use performed by the Depot/Seller or at a workshop authorized by Dentsply.
2. Before putting the Device into service, the directions of the Instructions for Use for proper use of the Device shall be followed in their entirety.
3. The special packaging shall be kept in a safe place, and in the event of any repairs, the Device shall be packaged as per the packaging instructions.

Section 5 Assignment and statute of limitations

1. If the Device covered by this Guarantee is sold to a third party, the terms and conditions of this Guarantee shall apply to the subsequent purchaser along with ownership of the Device.
2. Any claims under this Guarantee shall be subject to a limitation period of six (6) months ending no later than six (6) months following the end of the guarantee term.

Section 6 Neglected obligations

If the Buyer/Beneficiary fails to meet his or her obligations under Section 4 para. 1 and 2, the Manufacturer/Guarantor shall be free from all obligations under this Guarantee.

Section 7 Statutory warranty

This Guarantee does not cover the statutory warranty. The terms of the statutory warranty shall apply in addition to the term of this Guarantee. If the Buyer/Beneficiary fails to meet his or her obligations under Section 4 para. 1 and 2, the warranty shall be void, unless the Buyer can show that his or her failure to meet his or her obligations did not result in any damage.

Section 8 Governing Law

This Guarantee shall be covered by German law. The place of jurisdiction shall be Hanau am Main, Germany.

Notice: The Terms and Conditions of the device manufacturers shall apply.

11 Declaration of Conformity

EU Declaration of Conformity, EU Machines Directive 98/37/EG

DeguDent GmbH

Rodenbacher Chaussee 4
D-63457 Hanau-Wolfgang

We hereby declare that the :

Multimat®2 Touch
Multimat®2 Touch+Press

conforms with the basic requirements of the following EU directives:

- 1. EU Machinery Directive 98/37/EG**
- 2. EEC Directive on electrical equipment designed for use within certain voltage limits 73/23/EEG**
- 3. EEC Directive relating to electromagnetic compatibility 89/336/EEG**

The following harmonized standards were employed:

DIN EN 50081-1: 1993-03

DIN EN ISO 12100-1: 2004-04

DIN EN 50082-2: 1996-02

DIN EN ISO 12100-2: 2004-04


DIN EN 61010: 1994-03

EN 60204-1: 1998

Hanau, August 5, 2004

Dr. Udo Schusser
Research & Development

Torsten Schwafert
Director business area dental laboratory products

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Dentsply International
P.O. Box 872
570 West College Avenue
York, PA 17405
USA

Dentsply DeTrey S.A.
17, Michael Faraday
F-78180 Montigny Le Bretonneux
France

Dentsply Canada Ltd.
161 Vinyl Court
Woodbridge, Ontario L4L 4A3
Canada

Dentsply Argentina S.A.C.I.
Gral. Enrique Martínez 657/661
RA-1426 Buenos Aires
Argentina

Dentsply Australia PTY Ltd.
204-206 Gipps Street
Abbotsford, Victoria 3067
Australia

Dentsply Asia
23/F Gee Chang Hong Centre
65 Wong Chuk Hang Road
Aberdeen
Hong Kong

Dentsply Sankin K.K.
14-9, Yushima 3-chome
Bunkyo-ku, Tokyo 113 - 0034
Japan

Dentsply S.A. (Pty) Ltd.
197 Biggersberg Road
Glenvista Ext. 6
2058 Johannesburg
South Africa

Dentsply Industria e Comercio Ltda.
Rua Alice Hervé 86
Bengen
25665-010 Petropolis, R.J.
Brazil

Dentsply Dental Co., Ltd.
H6 Block A
TEDA International Incubator Centre
No. 5 Avenue
Tianjin Economic Development
Area 300457
Tianjin
China

Dentsply India pvt. Ltd.
W-86, Lower Ground
Greater Kailash-II
New Delhi 110048
India

Dentsply Korea
Hyun Woo Building 1st Floor
459-5, Dogok-dong
Kangnam-ku, Seoul
Korea

Dentsply México, S.A. de C.V.
Calzada Vallejo No. 846
Col. Industrial Vallejo
02300 México, D.F.

Dentsply Philippines Inc.
3/F Trans-Phil House
1177 Chino Roces
Avenue Car
Bagtikan Street
1203 Makati City
Philippines

Dentsply Russia Ltd.
Kozhevnikeskaya Street 10/2
6th Floor
113114 Moscow
Russia

Dentsply Ukraine
Krasnoarmeyskaya Street 912
Office 40
Kiev 01004
Ukraine
Tel. Kiev 4952226

Dentsply Taiwan
7F-11, 351,
Sec.2 Chung Shan Road
Chung-Ho, Taipei County
Taiwan R.O.C.

Dentsply Thailand Ltd.
23rd Floor Panjathani Tower
127/28 Ratchdapisek Road
Chongnonsee, Yannawa
Bangkok 10120
Thailand

Dentsply Vietnam
4th Floor,
6 Phung Khac Khoan Str.1
Ho Chi Minh City
Vietnam

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Dentsply DeTrey GmbH